



Magnusson entrepreneurial experience

- Co-founder and Chief Technical Officer, Resonant Optics Incorporated
- Co-founder and Chief Technical Officer, Resonant Sensors Incorporated
- Founder, Tiwaz Technologies LLC

RSI: provides next-generation optical sensor systems for pharmaceutical and biotech customers including COVID-19.



ResoSens™ bioassay system




Bionetic
microarray plates



Set Your **Label Free**

Commercial system



ResoSens™ system including ResoVu™ analysis software

- **Label-free** operation in a microarray format
- Automation compatible
- Windows 8 based, exportable CSV data files
- User selected plate layout and beam size—easy to use!


96- & 384-well microarray plates

- Bare and pre-activated
- Custom assay kits and/or reagents
- Method development and/or testing services
- Detect an array of analytes in a single well
 - custom multiplexed kits (4- or 8-plex)

<http://resonantsensors.com/>



96 wells
Scan <20 sec




384 wells


Bionetic

High-throughput ResoSens COVID-19 tests


- Patient **serum/blood/saliva** tests for antibodies against COVID-19 antigen
- Tests are performed using COVID-19 receptor binding domain (RBD) antigen.





Scientists exploring how coronaviruses like COVID-19 infect human cells have shown that the **SARS-CoV-2 spike (S) glycoprotein binds to the cell membrane protein** angiotensin-converting enzyme 2 (ACE2) to enter human cells.

S-protein=Key


ACE2=Lock



Patient serum testing for antibodies against the COVID-19 antigen (RBD)



Add saliva or serum sample




**Results measured in quadruplicate and averaged.
Reference buffer subtracted from data.**

**4 wells/sample; 24 tests in a single run; ~30 minutes/run;
<2 min/test**

WORKS GREAT!



New system design and color scheme



HEX color
0a014f
for metal
housing on
system and
incubator.

New project: LWIR optical components

Multi-layer film stacks versus single-layer resonance devices




Fig. 2.9. (a) An example of a Bragg stack. Thicknesses, $d(\text{SiO}_2) = 303.82 \text{ nm}$, $d(\text{TiO}_2) = 180.04 \text{ nm}$; refractive indices $n_c = 1.00$, $n(\text{SiO}_2) = 1.44$, $n(\text{TiO}_2) = 2.43$, $n_s = 3.46$; incident angle $\theta_{in} = 0^\circ$, (b) Calculated spectral response of a Bragg stack for TE polarized waves.




Fig. 2.10. (a) Configuration of a GMR high reflector. Thickness, $d = 490 \text{ nm}$; refractive indices $n_c = 1.00$, $n_H = 3.48$, $n_L = 1.00$, $n_s = 1.48$; grating period $\Lambda = 766 \text{ nm}$; filling factor $f = 0.726$; incident angle $\theta_{in} = 0^\circ$, (b) Calculated spectral response of a GMR device for TM polarized waves [from reference 14].

Beyond thin-film optics:


Resonant-grating-based optical component technology

NSF-SBIR Phase I 2023-2024

Optical component business ~\$5-10B/year

This year: Phase I is 275k, Phase II 1M at NSF=National Science Foundation

Table I: Thin-film multilayer reflector band coverage provided by key companies



The electromagnetic band from ~8 to 12 μm covers a region of atmospheric transparency important for long-range terrestrial imaging, spectroscopic applications, night-vision systems, and medical and industrial laser technologies. To utilize this spectral band for applications, effective components including reflectors, filters, and polarizers must be available.

Preliminary results LWIR optical components

